

LEVEL 1A

OFF-SITE CERTIFICATION & TRAINING



To prepare for an off-site class, please review and check off the items listed below. If your lab does not carry or cannot obtain the material or equipment listed, you will not be able to schedule an off-site class. Once your checklist is complete, please have a supervisor sign and date to acknowledge that all requirements have been met. This checklist will need to be submitted to the HMAC prior to any scheduling of an off-site class.

CONNECTION

- TV or projector for presentations (must be able to cast from Instructor laptop)
- Stable Wi-Fi connection

SETUP

- Clean lab with minimum interruptions and space for 6 people
- Computer or laptop for CERTIFICATION CLASSES ONLY (1/person)
- Testing or lecture room with plenty of seating space for 6 people

EQUIPMENT

- Tex-207-F** Texas Gyratory
- 960-gram pre-heated samples (2/person)

- Tex-207-F** Scale
- Water Bath Tank with Heater Circulator
- Suspension Basket
- Calipers
- Damp towel
- (1) Laboratory Molded Sample
- (1) Road Core
- Metal Pycnometer
- Vibrating Table
- Vacuum Pump
- *(2) 1200-gram samples

- Tex-241-F** Superpave Press
- 2400-gram pre-heated sample (2/person)

- Tex-236-F** Ignition Oven
- Basket
- *(1) 1200-gram sample to burn sample
- (1) after-burnt sample

- Tex-200-F** Burnt Sample from Tex-236-F
- Sieve Scale
- Sieve Pans
- Tools

**Use 3X amount for training*

Supervisor Name (print)

Supervisor Signature

Date